

KP Series Cannulas

for Organ Perfusion

Designed for Transplant Procedures

- FDA Cleared
- French-sized
- Color-coded
- Suture bump for secure suturing
- Adjustable suture bump (16F and up) for additional secure fixation
- Flexible kink-resistant connector
- High flow
- Rubber connectors are firm to grip but flexible
- Back table adaptor designed for removal of perfusion fluid residues in organs

To order or for more information contact your local representative or orders@B2LL.com

French-Sized Cannulas and Tubing



KP10-16



KP12-19



KP16-22-26



KP20-28-32



KP25-34-40

Cannula Product Information

ORDER#	SIZE (French)	LENGTH (mm)	SUTURE BUMP	COLOR	PER CASE
KP10-16	10	305	—	Black	10
KP12-19	12	305	—	Clear	10
KP16-22-26	16	305	Yes	Orange	10
KP20-28-32	20	305	Yes	Yellow	10
KP25-34-40	25	305	Yes	Blue	10
KP25-34-K	25 (short)	115	—	Orange	10
KP25-34-K-B	25F (short-bent)	95.5	—	Blue	10

Accessory Product Information

ORDER#	DESCRIPTION	PER CASE
BTL-173	Quad-Spiked Organ Perfusion Set	15
BTL-174	Dual-Spiked Organ Perfusion Set	15
BTL-177	Portable Tubing	25

Now Available



115 mm cannulae KP25-34-K: ~60 mm



25F Short-bent KP25-34-K-B



Back Table Adaptor KP16-32A: ~130 mm